

Part Number: BE-35-04-110	Date: 27OCT20	Proofread No.: P04
Designer: RK	Page: 1	
Colour Reference: <span style="background-color: red; color: white; padding: 2px;">P185U</span>		

**AVDB3WH1E**  
**Baxter**

ฉลากถุงบรรจุยา ตำรับยาชื่อ OLIMEL N12E

ขนาดบรรจุ 1000 ml

**1,000 mL**

Bag

Nitrogen g	12.0
Total kcal	950
Non-protein calories	640

**Central Vein Only**

**With Electrolytes**

“百特” 歐立美N12E 輸注乳液  
**OLIMEL N12E**

**12 g/l nitrogen 950 kcal/l with electrolytes Emulsion for Infusion**

僅供中央靜脈輸注

**For central intravenous administration**  
**1,000 mL of the final mixture provides the following:**

Total calories approx.	950 kcal
Non-protein calories	640 kcal
Nitrogen	12.0 g
Non-protein calories / nitrogen ratio	53 kcal/g
Glucose	73.3 g
Lipids	35.0 g
Glucose / lipid calories ratio	45/55
Amino acids	75.9 g
Glucose calories	290 kcal
Lipid calories	350 kcal
Lipid / total calories	37%
Electrolytes:	
Sodium	35.0 mmol
Potassium	30.0 mmol
Magnesium	4.0 mmol
Calcium	3.5 mmol
Phosphate	15.0 mmol
Acetate	70 mmol
Chloride	45 mmol
pH	6.4
Osmolarity	1,270 mOsm/L

**Formula for 1,000 mL after mixing:**

Refined olive oil + refined soybean oil 35.00 g, Alanine 10.99 g, Arginine 7.44 g, Aspartic acid 2.20 g, Glutamic acid 3.79 g, Glycine 5.26 g, Histidine 4.53 g, Isoleucine 3.79 g, Leucine 5.26 g, Lysine (equivalent to lysine acetate) 5.97 g (8.43 g), Methionine 3.79 g, Phenylalanine 5.26 g, Proline 4.53 g, Serine 3.00 g, Threonine 3.79 g, Tryptophan 1.26 g, Tyrosine 0.20 g, Valine 4.86 g, Sodium acetate, trihydrate 1.5 g, Sodium glycerophosphate, hydrated 3.67 g, Potassium chloride 2.24 g, Magnesium chloride, hexahydrate 0.81 g, Calcium chloride, dihydrate 0.52 g, Glucose (equivalent to glucose monohydrate) 73.33 g (80.67 g)

**Excipients:**

purified egg phosphatide, glycerol, sodium oleate, sodium hydroxide, water for injections, glacial acetic acid, hydrochloric acid.

After reconstitution: Emulsion for Infusion 1,000 mL in a 3-compartment bag

Glucose solution 27.5% : 267 mL  
Amino acid solution 14.2% : 533 mL  
Lipid emulsion 17.5% : 200 mL

百特醫療產品股份有限公司  
台北市敦化南路二段95號28樓  
製造廠: Baxter S.A

廠址: Boulevard René Branquart, 80, 7860 Lessines, Belgium

適應症、用法用量及注意事項:  
請詳見仿單  
請勿冷凍

Manufactured by: Baxter SA, 7860 Lessines, Belgium  
Imported by: Baxter Healthcare (Thailand) Co, Ltd., Bangkok, Thailand  
Thailand Reg. No. 2C xx/6x

For single use only

Read the package leaflet before use

Keep out of the sight and reach of children

Do not freeze

Store in the overpouch

Do not infuse before, simultaneously or after administration of blood through the same equipment

Use only if the amino acid and glucose solutions are clear, the lipid emulsion and the final mixture after reconstitution are homogenous and with a milky appearance and the container undamaged

Do not connect in series

Do not infuse without opening the seals and mixing the contents of the 3 compartments

After reconstitution, immediately administer the product

**500**

ยาอันตราย



Mfg. Date:

**AVDB3WH1E**

(01)00085412676999



ยาสิ้นอายุ  
Exp:

ขอรับรองว่าข้อความบนฉลากของ  
ภาชนะบรรจุทุกชนิดและหีบห่อ  
ถูกต้องตรงตามมาตรฐาน ๒๕ (๓)  
หรือ ๒๗ (๒) และ (๓) แล้วแต่กรณี

*John S.M.*

**POM**

BE-35-04-110

Lot:

