

ฉลากกล่องยา 700 mg

**INVEGA HAFYERA®**  
prolonged release suspension  
for intramuscular injection

**700mg\***

SHIP AND STORE THIS SIDE UP LAY FLAT

Administer every 6 months

Shake syringe with the syringe tip cap pointing up VERY FAST for at least 15 seconds, rest briefly, then shake again for 15 seconds

Do not store above 30°C. Keep out of reach of children.

Single-dose prefilled syringe. Use entire contents of syringe.

FOR GLUTEAL INTRAMUSCULAR INJECTION ONLY

CONTENTS: 1 single-dose prefilled syringe and 1 needle (a 20G,1.5-inch thin wall safety needle).

For directions for use, refer to package insert.

Keep out of the sight and reach of children.

Do not store above 30°C.

Manufactured by Janssen Pharmaceutica NV, Beerse, Belgium

Imported by Janssen-Cilag Ltd., Bangkok, Thailand

Please see full prescribing information and company contact information at <https://ind.fda.moph.go.th>

Reg. No. TC 15036/66(N)

สำหรับใช้ฉีดเข้ากล้ามเนื้อ  
ใช้เข็มที่แถมมา  
อย่าให้เด็กเห็น

© J-C 2024

**INVEGA HAFYERA®**  
prolonged release suspension  
for intramuscular injection

**700mg\***

SHIP AND STORE THIS SIDE UP LAY FLAT

Administer every 6 months

Shake syringe with the syringe tip cap pointing up VERY FAST for at least 15 seconds, rest briefly, then shake again for 15 seconds

Do not store above 30°C. Keep out of reach of children.

Single-dose prefilled syringe. Use entire contents of syringe.

FOR GLUTEAL INTRAMUSCULAR INJECTION ONLY

Each injection must be administered only by a healthcare professional.

Shake before using.

Do not store above 30°C. Keep out of reach of children. Rx Only

ดูก่อนใช้เพื่อความปลอดภัยกับยาของคุณมากที่สุด โดยใช้แอปพลิเคชัน Scan Matrix หรือ <https://epi.jnj>

Scan for latest E-leaflet by using Scan Matrix application or <https://epi.jnj>

8 850583 003330

**Johnson & Johnson**

ขอรับรองว่าข้อมูลความมั่นคงปลอดภัยของเอกสารบรรจุภัณฑ์นี้ดูแลและพิมพ์โดย Chayapitcha Saesue (ช) และ (ค) แล้วแต่กรณี

Digitally signed by Chayapitcha Saesue  
DN: c=US, o=JNJ, ou=Subscribers, email=chayapitcha.saesue@us.jnj.com, cn=Chayapitcha Saesue  
Reason for latest to the accuracy and integrity of this document:  
Date: 2024.07.24 09:54:05 +07'00'  
Unique identifier:  
20240713.20064

AUTOMATED LEGEND (CC-PO)		Product Artwork	
Box BYANNLI	PCC: MU_12352573	Latest Modification Date: July 11, 2024 04:58 PM	Version: 1
Template Ref. Name: BeerselBox953S5	Template Dimensions: 124 x 39 x 189	Artwork Operator: V_PD	
Tech. Info/Spec:		FINIS 915	<input type="checkbox"/> Decar (not to be printed)
Market: TH	dummy	FINIS 803	<input type="checkbox"/> Wash free
Matr. ID Code: MU_12352573.ai		FINIS 2947	<input checked="" type="checkbox"/> Wash free & Used for glass
Replace PCC:		Black	Perforation: .....
Font Type: Johnson Text, Brouallia New			